# MD Argon

Argon-Druckminderer für MD Argon Geräte / Argon Pressure Reducer for MD Argon Units / Réducteur de pression pour MD Argon / Reductor de presión para MD Argon / Riduttore di pressione per MD Argon

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Seite 3

DEUTSCH

## Operating Manual

Page 13

ENGLISH

## Mode d'emploi

Page 23

FRANCAIS

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Página 33

ESPAÑOL

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Página 43

ITALIANO

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5.3.1 Opening ..... 18

5.3.2 Closing..... 18

5.3.3 Adjusting the mechanical style pressure switch adjustment screw ..... 18

**6 Cleaning..... 18**

**7 Maintenance..... 19**

7.1 Changing gasket ring..... 19

7.2 Changing filter screw..... 19

**8 Storage ..... 19**

**9 Technical data ..... 19**

9.1 Pressure reducer..... 19

9.2 Mechanical style pressure switch..... 20

**10 Servicing and manufacturer address..... 20**

# 1 Symbol Descriptions

Please read these instructions carefully and keep them for future reference. Improper handling can cause personal injury or damage the product.

Symbol	Definition
	Batch Code
	Medical Device
	Catalog number
	Unique Device Identifier
	GS1 standard
	Caution!
	Non-sterile
	Consult Instructions for Use
	CE-marking
	Manufacturer
	Production date
	Keep dry
	Keep away from sunlight

## 2 Cautions

The following safety precautions must be observed:

Incorrect use of the pressure reducer can cause serious damage.

Run the pressure reducer only at the intended pressure and with the intended gas. The use of improper gas or at improper pressure might cause risks of explosion.

Protect the pressure reducer from dust, oil and other sources of dirt.

Do not use the pressure reducer if it is not in perfect working condition.

When you draw gas, the cylinder must be placed upright and protected from falling.

Danger of explosion. At startup and also during use, hands and tools must be free from oil and grease.

## 3 Intended use

This product delivers Argon gas to **MDV touch 300 units**. It reduces the gas pressure of the source at the required level and maintains it.

## 4 Product description

- High pressure gauge (1)
- Mechanical style pressure switch (2)
- Pressure adjusting screw (3)
- Relief valve (4)
- Outlet connection (5)
- Body (6)
- Gasket ring (7)
- Coupling nut (8)
- Inlet connection (9)
- Filter screw (10)
- Adjusting screw mechanical style pressure switch (11)


Figure 1

Figure 2

## 5 Commissioning and use

### 5.1 Connection to gas cylinder

The pressure reducer is constructed for gas cylinder connections in line with DIN 477.

- At every gas cylinder change, initially check that the gasket ring (7) is intact.
- Before connecting the pressure reducer to the gas cylinder, make sure the openings of the gas cylinder and the pressure reducer are free from dirt. Otherwise blow out with compressed air.
- Turn the pressure adjusting screw (3) anti-clockwise as far as it will go. Do the same with the valve of the gas cylinder.

- Put the inlet connection (9) of the pressure reducer on the opening of the gas cylinder. Align the high pressure gauge (1) upright and close the coupling nut (8) manually. If necessary again align the high pressure gauge.
- Tighten the coupling nut by using a 30 mm flat wrench.

Stay away from heat sources and easily inflammable materials during the commissioning.

Do not direct the opening of the gas cylinder on yourself or other persons.

## 5.2 Connection to HF-Device

### Ports of the HF-Device


Figure 3

- gas pressure measure line 2 (12)
- gas connection 2 (13)
- gas pressure measure line 1 (14)
- gas connection 1 (15)

- Plug in the mechanical style pressure switch (2) leaving cable in gas measure line 1 or 2 (14 or 12).
- Turn the lock clockwise until it is engaged. You can hear that at the click-sound.
- Plug in the outlet connection (5) leaving tube with its adapter at the end to gas connection 1 or 2 (15 or 13). The lock engages by itself (Figure 4 and 5).


Figure 4

Figure 5

## 5.3 Use

### 5.3.1 Opening

- Slowly open the gas cylinder valve. The high pressure gauge (1) will show you the cylinder pressure.
- Slowly open the pressure adjusting screw (3) as far as it will go.

**Opening the gas cylinder valve too quickly might cause the gauges to malfunction.**

### 5.3.2 Closing

- Close the gas cylinder valve.
- Release the gas out of the pressure reducer.
- Turn the pressure adjusting screw (3) anti-clockwise till it is completely closed.

### 5.3.3 Adjusting the mechanical style pressure switch adjustment screw

The mechanical style pressure switch fulfills the purpose to provide the control system of the HF-device with information, whether the required pressure is on or not. For that you can set a minimum pressure with the aid of the adjustment screw. Pay attention that the minimum pressure is always set under 4 bar, otherwise the gas supply will be inhibited.

To set the minimum pressure, perform the following instructions:

- Turn the adjustment screw (11) clockwise with the aid of a slotted screwdriver, to increase the pressure.
- Turn the adjustment screw (11) anti clockwise with the aid of a slotted screwdriver, to decrease the pressure.
- The adjusted pressure can be seen on the scale.

## 6 Cleaning

In case you need to clean the pressure reducer use a dry, dust-free cloth to wipe the item. For stubborn soiling use a damp cloth.

**Avoid water inside the pressure reducer. It might cause risks of explosion!**

## 7 Maintenance

Hands and tools must be free from oil and grease during maintenance

### 7.1 Changing gasket ring

At a high level of abrasion perform the following instructions:

- Lever the gasket ring with the aid of a slotted screwdriver and take it off.
- Pull the fresh gasket ring into place without any aid of tools.

### 7.2 Changing filter screw

At a high level of soiling, but at least every two years, perform the following instructions:

- Release the filter screw with the aid of a 6 mm Allen wrench.
- Hold the pressure reducer so that the inlet connection (9) is facing to the ground.
- This prevents dirt getting inside.
- Manually turn out the filter screw.
- Manually turn in the fresh filter screw.
- Tighten the filter screw with the aid of a 6 mm Allen wrench.

## 8 Storage

The pressure reducer must be stored dry and dust-free. The surrounding air must be non-condensing. It is recommended to keep the wrapping, so the item can be storage in it when required.

## 9 Technical data

### 9.1 Pressure reducer

Inlet pressure P1 [bar]	max. 230
Outlet pressure P2 [bar]	4
Delivery flow Q1 [m <sup>3</sup> /h]	< 1
Operating temperature [°C]	-20 bis +60

## 9.2 Mechanical style pressure switch

Model	IS 1000		
Proof pressure [bar]	10		
Max. pressure [bar]	7		
Operating pressure range [bar]	1 bis 6		
Hysteresis [bar]	0,8		
Ambient and gas temperature [°C]	-5 bis +60		
Max. contact capacity	2 VAC AC / 2 W DC		
Voltage	24 V AC/DC	48 V AC/DC	100 V AC/DC
Max. operating current	50 mA	40 mA	20 mA

## 10 Servicing and manufacturer address

If you require a paper copy of these instructions for use, please contact Micromed Customer Service using the contact details provided. The paper copy of the instructions for use will be provided to you within seven calendar days of receipt of your request. Alternatively, you can use the attached QR code.

MICROMED Medizintechnik GmbH  
 Eisenbahnstraße 84  
 78573 Wurmlingen Germany  
 Tel.: 07461 / 96 48 55-0  
 Fax: 07461 / 96 48 55-95  
 E-Mail: [info@micromed.com](mailto:info@micromed.com)  
<https://www.micromed.com/downloads/>

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Eisenbahnstr. 84  
78532 Wurmlingen / Germany  
Tel. +49 (0) 7461 96 48 55-0  
Fax +49 (0) 7461 96 48 55-95  
info@micromed.com  
www.micromed.com

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