MDS Smoke Evacuator

Rauchgasabsaugsystem / Smoke evacuation system

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Seite 3

Operating Manual

Page 22

DEUTSCH

ENGLISH

Inhalt

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Contents

1	Symbol descriptions	23
2	Introduction	24
3	Scope of delivery	24
3.1	Optional accessories	25
4	Scope	25
4.1	Inteded Use	25
4.2	Contraindications	25
5	Structure and basic components of the system	25
5.1	Main menu	26
5.2	Selection of operation mode	27
5.3	Setting the time for the delayed motor stop	27
5.4	Selection of suction activation method	28
5.5	Additional Symbols on the MDS Smoke Evacuator	29
5.6	Filter life indication	29
5.7	Filter replacement	29
6	Specifications	30
7	Accessories	31
7.1	Liquid separators for flue gas extraction units	31
7.1.1	Construction and connections	32
7.1.2	Intended use	32
7.1.3	Scope	32
7.1.4	Principle of operation	32
7.1.5	Commissioning	33
7.1.6	Operation	33
7.1.7	Exchange	33
7.1.8	Safety instructions	33
7.1.9	Storage	33
8	Safety guidelines	33
9	Servicing and manufacturer address	35
10	Cleaning and disinfection	35
11	Maintenance and service	36
11.1	Repair	36
11.2	Replacement of fuses	36
11.3	Technical safety check	36
11.4	Transport	36
11.5	Environmental directives	37

1 Symbol descriptions

Symbol	Definition
	CE-labelling
	Attention
	Manufacturer
	Reference code
	Serial number
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
	Instrucions for use
	WEEE-labelling
	Green point

2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

The MDS Smoke Evacuator is designed to remove and filter smoke and aerosols from the surgical field generated during electrosurgical and laser procedures.

The system removes surgical smoke generated in the operating field using a powerful vacuum motor. The smoke is aspirated from the surgical site through a suction hose and directed into a multi-stage filtration system, where it is effectively treated and purified.

The filter used is designed as a single-use, fully encapsulated filter. This design allows for easy, quick, and safe installation and replacement of the filter. The complete encapsulation reliably protects medical personnel from potential contamination during filter replacement.

The air filtration process consists of four coordinated filtration stages:

1. Pre-filter
In the first stage, coarse particles and larger aerosols are captured and removed, thereby reducing the load on the subsequent filtration stages.
2. ULPA filter (Ultra Low Penetration Air)
The second stage consists of a high-efficiency ULPA filter that captures particles and microorganisms in the size range of 0.1 to 0.2 micrometers with an efficiency of 99.999%.
3. Activated carbon filter
In the third stage, activated carbon is used to remove gaseous pollutants, odors, and volatile organic compounds from the surgical smoke.
4. Glass fiber post-filter
The fourth stage consists of a filter medium made of woven glass fiber fabric. This stage minimizes the passage of fine activated carbon particles, further enhancing the cleanliness of the exhausted air.

The MDS smoke evacuator is equipped with a modern touchscreen display that enables intuitive operation. The LCD screen allows the user to select and monitor the programming mode, motor power, airflow rate, activation method, as well as the filter status and remaining filter life.

The device can be activated flexibly via the on-screen control, a pneumatic foot pedal or an automatic RF-sensor that enables demand-controlled operation.

3 Scope of delivery

155-100-001	MDS Smoke Evacuator
155-102-001	Filter for MDS Smoke Evacuator - 3-Port System
151-002-200	Power cable
	Instructions for Use

3.1 Optional accessories

108-241-324	Disposable electrode handle with 2 buttons, with spatula, 3 m cable, with smoke extraction Ø 2.4 mm, sterile packaging
152-101-001	Pneumatic foot switch
155-102-001	Replacement filter
155-103-001	HF sensor connection cable
155-104-001	Liquid separator

4 Scope

The product is used in electrosurgery. It may only be used by trained and qualified professionals.


4.1 Inteded Use

A device for extracting smoke and aerosols. The device is used during electrosurgical procedures to extract and filter smoke, thereby ensuring visibility.

4.2 Contraindications

No contraindications are known.

5 Structure and basic components of the system

Front view of the MDS Smoke Evacuator

Located on the front of the unit are:

- (1) Filter system unit
- (2) HF-sensor input
- (3) HMI touch display
- (4) Locking mechanism for the filter system unit
- (5) AIR footswitch input
- (6) On/Off main switch


Rear view of the MDS Smoke Evacuator

Located on the back of the housing are:

- (7) Power input module with build-in fuses holder
- (8) Casing grounding
- (9) Label

5.1 Main menu


Main menu of MDS Smoke Evacuator

- (10) Filter life indication. See Chapter 5.6
- (11) Indication/Button for selecting the operating mode (See Chapter 5.2).
- (12) Indication of the selected power of the motor
- (13) Button to increase the power of the motor, in steps from 1 to 10, each step increases the power by 10%
- (14) Motor power reduction button, in steps from 10 to 1, each step reduces power by 10%
- (15) Suction enable/disable button. When the turbine is not running, the button is white
- (16) Indication of the state of the turbine - animation at the start of suction
- (17) Button for setting the delayed time for turning off the turbine (See Chapter 5.3)
- (18) Indication/Button for selecting the activation method
- (19) Indication of the level of airflow through the filter

5.2 Selection of operation mode

To access the "Operation Mode Selection" menu, simply press the button "11".


- **(20)** "Laparoscopic procedure" mode - if the system is connected to a surgical cannula during a laparoscopic operation
- **(21)** "Electrosurgical Mode" - if an electrocautery handle is used or an integrated smoke evacuation handle is used
- **(22)** "Open procedure" mode - an open procedure in which a smoke evacuation hose is used
- **(23)** Exit button - This button allows you to exit to the main menu and save the selection

5.3 Setting the time for the delayed motor stop

The system has a feature that enables users to set a time for turning off the motor after deactivating the HF-energy. This allows the motor to continue running for a specific time to clear out any residual smoke from the operative area, after switching off by either the onscreen button, the footswitch or the HF-sensor.

To access the settings menu, press button 17 from the main menu. This will allow you to set the time between 0 and 10 seconds.

5.4 Selection of suction activation method

Selecting and activating the suction method

Press button 18 in the main menu to select the desired suction method. There are two ways to activate the suction:

- Activation via a pneumatic foot switch
- Activation via an HF sensor

Suction is generally activated using button 15 in the main menu, provided that the selected method allows this.

Activating suction with a pneumatic foot switch

When activating via a pneumatic foot switch, the connection hose is plugged into socket 5 (connection port for the pneumatic foot switch) on the front panel.

To do this, select the “Pedal start” function in the main menu.

Suction is switched on and off by pressing the pneumatic foot switch.

Suction activation with HF-sensor

When activating via an HF sensor, the connection cable is plugged into socket 2 (HF-sensor input) on the front panel.


To do this, select the “
function in the main menu.

The HF-sensor is clipped onto the connection cable of the electrode handpiece.

Suction starts automatically as soon as the HF-sensor detects HF energy on the connection cable and stops automatically as soon as no HF current is flowing to the electrode handpiece.

Note:

When HF-sensor control is activated, neither the pneumatic foot switch nor button 15 can be used to control the suction.

The motor will be operated by pedal

The motor will be controlled by the HF-sensor

5.5 Additional Symbols on the MDS Smoke Evacuator

Symbol	Meaning
Symbols on the frontplate	
	Consult instruction for use
	Footswitch plug-in connector
	HF-sensor plug-in socket
Symbols on the backplate	
	Grounding

5.6 Filter life indication

Symbol	Definition
	New filter
	The filter has been working for more than 16 hours.
	The filter has been working for more than 33 hours.
	The filter is not inserted or its life of 40h has expired

5.7 Filter replacement

The filter system unit is designed as a single-use filter and has a maximum operating life of 40 operating hours. After this period, the filter performance is exhausted and safe operation of the smoke evacuation system is no longer possible. The filter system unit must therefore be replaced.

Before replacing the filter, ensure that the smoke evacuation system is switched off and disconnected from the power supply. This prevents potential injury and unintended aspiration of air or particles during the filter replacement process.

To remove the used filter system unit, first release the locking mechanism of the filter system unit. After unlocking, the filter unit can be pulled forward and removed from the smoke evacuation system. Due to the complete encapsulation of the filter unit, medical personnel are protected from potential contamination during this procedure.

The new filter system unit is then inserted from the front into the smoke evacuation system. Ensure that the filter unit is correctly aligned and fully inserted into the device without tilting or jamming. After insertion, the filter system unit must audibly and tactilely engage in the designated locking mechanism.

Finally, verify that the filter unit is firmly and securely seated in the smoke evacuation system. Operation of the device is permitted only when the filter system unit is properly locked in place. Only after this verification may the smoke evacuation system be reconnected to the power supply and put back into operation.

The removed filter system unit must be disposed of in accordance with applicable clinical and legal regulations for contaminated medical disposable products.

Filter unit

6 Specifications

Maximum flow setting in "Open procedure" program		
7/8" (22mm)	35 CFM	990 LPM
3/8" (9,5mm)	4.5 CFM	130 LPM
1/4" (6,4mm)	2 CFM	57 LPM
Device Dimensions (H x W x D)	370mm x 144mm x 375mm	
Weight	10kg	
Filter Dimensions (H x W x D)	128mm x 112mm x 243mm	
Noise Level, dBA	55dB	
Footswitch Activation	Standard pneumatic	
Remote Control Activation	RF sensor in range 350 - 600kHz	
Flow Control	Variable 1 - 10 Steps, in step of 1	
OFF Time Delay Control	Setting from OFF to 10sec, in step of 1 sec	
Filter Status		

Green	New
Orange	Middle, >16h
Red	High, >33h
Violet	Service, 40h
Vacuum level	Yes, bargraph
Choice Of Three Programs	
Open surgical Tubing	100% of Maximal flow
Electrosurgical Pencil Attachment/Integrated Smoke Evacuation Pencil	85% of Maximal flow
Laparoscopic Tubing	75% of Maximal flow
Choice of activation method	Footswitch, Remote control
Power supply	110 – 260 Vac, 50 – 60 Hz
Input power	Max 500VA (PF = 0.5) Standby 6VA (PF = 0.4)
Fuses	2 x T6.3A
Construction	EN 60601-1
Protection class	I
Grounding terminal	Yes
Motor static suction	21 kPa


7 Accessories

To avoid incompatibility and unsafe operation, use only the manufacturer's recommended accessories for working with the MDS Smoke Evacuator.

Follow the instructions to clean and prepare the accessories.
Before each sterilization, check the connecting cables and handles for damage. Check for broken insulation, especially the electrode cables and the accessories used for endoscopic manipulations.


7.1 Liquid separators for flue gas extraction units

7.1.1 Construction and connections

There are three hose connections of different diameters on the front of the liquid separator. These are used to connect an electrode handle with an integrated flue gas extraction system and allow the use of electrode handles with different hose diameters. When in use, care must be taken to ensure that the unused hose connections are securely closed with the cap.

On the rear of the liquid separator there is a connection socket through which the liquid separator is inserted into the flue gas extraction unit. The liquid separator is secured in the flue gas extraction unit by both form and force fit and is securely positioned during operation.

7.1.2 Intended use

The liquid separator is designed to separate liquids and condensates from the smoke gas extracted during electrosurgical procedures. It protects the smoke evacuation unit and its filter system from the ingress of liquids and contributes to safe operation. The liquid separator is intended for single use and must be disposed of after use.

7.1.3 Scope

The liquid separator is intended exclusively for use with MDS flue gas extraction units designed for this purpose and approved by the manufacturer.

7.1.4 Principle of operation

The extracted flue gas is passed through the liquid separator. The extracted liquids and condensates are retained in the collection tank by the design of the system. The pre-cleaned flue gas is then fed into the flue gas extraction unit's filter system.

7.1.5 Commissioning

1. Check the liquid separator for visible damage before use.
2. Fit the liquid separator between the extraction hose and the flue gas extraction unit, ensuring it is oriented in the correct flow direction.
3. Ensure that all connections are secure and leak-free.

7.1.6 Operation

- Always use the liquid separator whenever liquids or condensate are likely to be drawn in.
- Monitor the fill level of the collection tank during operation.
- Do not exceed the maximum permissible fill level.

⚠ Warning:

An overfilled liquid separator can impair extraction performance and cause damage to the equipment.

7.1.7 Exchange

- Disposable liquid separators must be disposed of after use in accordance with the applicable hygiene and waste disposal regulations.

7.1.8 Safety instructions

- Do not use the liquid separator if it shows any visible damage, such as cracks or breaks.
- Do not operate the flue gas extraction unit without the liquid separator correctly fitted, where one is required.
- Modifications, misuse or the use of unauthorised accessories are not permitted.

7.1.9 Storage

Store the liquid separator in a clean, dry place, protected from direct sunlight, in accordance with the manufacturer's instructions.

8 Safety guidelines

Micromed Medizintechnik mandates that the medical device should only be operated under the supervision of authorized or qualified personnel. The healthcare professional and surgeon are required to have adequate training and familiarity with the fundamental principles, application rules, and potential risks associated with electrosurgery to ensure the safety of patients, personnel, and equipment.

This system is intended to be used exclusively for medical purposes. For safety reasons, the connection to the mains must be made through properly protected sockets, using the manufacturer's main connecting wires and plugs, or equivalent quality components in accordance with VDE 0107. To ensure safety, extension cords and splitters should not be used.

The device comes with a ground potential balance cable socket, which must be connected to an appropriate connection in the operating room to avoid any static electricity when the apparatus is repositioned.

	Repair of the appliance should be carried out only by Micromed Medizintechnik or by a representative authorized by Micromed Medizintechnik.
	Appliance not to be used in areas with explosion risk.
	The device, incl. connection cables to be checked for defects before each use.
	Defects to be removed promptly.
	Unplug the unit from the electrical outlet before checking the system components.
	Make sure flammable materials used as cleaners, disinfectants, or tape removers have evaporated or been completely removed before starting the motor.
	The use of non-flammable detergents for cleaning and disinfection is recommended.
	Do not block either the pipe or the filter. If they become clogged or air movement through the filter is significantly restricted, the motor may overheat and cause damage to the unit.
	The smoke evacuator creates a strong vacuum. Adjust the airflow and the position of the inlet end of the tube to prevent injury to the patient and to prevent aspiration of surgical materials and surgical specimens.
	If the smoke evacuator is activated while the airflow is set to high speed, it may cause a sudden, strong suction. Check the airflow setting before starting work to prevent patient injury and prevent aspiration of surgical materials and surgical specimens.
	To maximize patient safety, tubing or suction tips should not come into direct contact with tissue. Failure to do so may result in injury to the patient.
	MDS Smoke Evacuator filters and disposable accessories are completely disposable. Please dispose of in accordance with local principles or regulations and hospital policy.
	Only filters provided by the manufacturer have been proven to be compatible with the MDS Smoke Evacuator. Do not use other filters with this system.
	Careful installation of hoses, adapters, and suction boxes is essential to prevent the motor from overheating and to avoid voiding the unit's warranty.
	It is imperative that the installation of this equipment is done correctly to prevent any blockage of the intake and exhaust openings located at the bottom of the system. Failure to do so may result in reduced performance, damage and/or render the system inoperable, which will void the warranty. Please ensure that the installation is carried out in accordance with the manufacturer's guidelines to avoid any such consequences.

	<p>Attention: Please note that this device is NOT suitable for liquid aspiration. In order to avoid any damage to the device and to ensure optimal performance, it is important to install liquid collection devices if there is a possibility of liquid coming into contact with the filter. Failure to install a liquid collection device may result in filter clogging and electrical damage.</p>
	<p>The filter should be changed according to the filter life guidelines. Filter life is determined by time only, not filter performance.</p>
	<p>Changes or modifications not expressly approved by Micromed Medizintechnik could void the user's authority to operate the equipment.</p>
	<p>If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.</p>
	<p>Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.</p>

9 Servicing and manufacturer address

Please contact your supplier in case of a problem that is insoluble for the operator or if the service is necessary. When installing the device, it informs you about the possibilities for technical maintenance and service at the specific point.

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

Micromed Medizintechnik GmbH
 Eisenbahnstraße 84
 78573 Wurmlingen / Germany
 Tel. +49 7461 96 48 55 - 0
 Fax. +49 7461 96 48 55 - 95
 E-Mail: info@micromed.com
 Web: www.micromed.com

10 Cleaning and disinfection

	<p>When cleaning the device: Before cleaning, disconnect the device from the mains. Before connection the appliance to the mains, make sure that all flammable cleaning agents have completely evaporated.</p>
---	---

Follow local hygiene guidelines.

- Clean all surfaces of the device (including the front) with alcohol-free cleaners and disinfectants.
- Wipe the LCD screen with a sponge or cloth soaked in clean water. Using aggressive cleaning agents will damage the display.

- Dry the device with a clean, lint-free cloth.
- Do not clean or disinfect surfaces of the case/front panel with flammable or explosive products. If their use is unavoidable, flammable or explosive products must be completely evaporated before switching on the device.
- The manufacturer recommends products with a pH between 9 and 10, e.g. MediClean forte by Dr. med. Weigert.
- Do not immerse the device in liquid.
- Do not sterilize the device with steam, ethyl oxide or other means.

11 Maintenance and service

11.1 Repair

Upon service: the repair of the device should be carried out only by Micromed Medizintechnik or by an authorized representative.

Do not open the device housing. The MDS Smoke Evacuator does not contain any user-replaceable parts. The system may only be repaired by a qualified service engineer, who is specially trained to repair the MDS Smoke Evacuator.

Interventions of other persons can lead to cancellation of the guarantee and pose a risk for health.

11.2 Replacement of fuses

The fuses are located on the back of the device in a socket (7) on the rear panel, where the mains cable is also connected.

The fuses must be replaced as follows:

- Use a small screwdriver to remove the fuse holder.
- Replace defective fuses with the type described in technical data table.
- Put the fuse holder into its slot.

11.3 Technical safety check

To ensure proper and safe operation of the device, it is recommended to be inspected every 12 months by Micromed Medizintechnik or a Micromed Medizintechnik-approved service representative. At least every 12 months an inspection should be carried out as follows:

- Confirming the availability of instructions and user manual.
- Inspect the device and accessories for mechanical damage.
- Confirming the legibility of safety labels.
- Confirmation of secure and correct attachment of the unit.
- Confirmation that devices that are not functionally and operationally safe will no longer be used.

11.4 Transport

Prevent mechanical damage and moisture.

If the unit has been stored for a long period in a cold environment, wait before removing the protective packaging until the unit has reached room temperature.

During transportation, the standard safety measures additionally apply.
Only ever transport the unit in its original packaging.

11.5 Environmental directives

From the implementation date of the European directive 2002/96/EU as national legislation, the following regulations apply:

- Electrical and electronic units must not be disposed of with normal household waste.
- The consumer is required by law to return electrical and electronic devices at the end of their service life to designated public collection points or to the place of purchase.

Disposal details are set out by the applicable local legislation in the relevant country. The symbol on the product, the instructions for use or the packaging indicates the relevant applicable legislation. By sending end-of-life devices for reuse, material recycling or other forms of utilization, you are making a major contribution to environmental protection.

Please also send the unit packaging for environmentally responsible recycling after the end of the product's life.

12 Warranty card

Warranty card

We grant a on all electrosurgical devices

2-year warranty

We grant a warranty in accordance with the statutory and country-specific regulations on all electrical medical products (verification required in the form of an invoice or delivery note). The warranty period is two years. In the event of material or production errors, we provide free repair of the unit. Any attempt to carry out independent repairs will render the warranty cover null and void.


Damage caused by incorrect handling, overloading or normal wear and tear is not covered under the terms of the warranty.

In the event of a warranty claim and resulting repair, please enclose a copy of the invoice or delivery note (or lot number) with the unit.

Manufacturer and aftersales service address:

Micromed Medizintechnik GmbH
Eisenbahnstr. 84
78573 Wurmlingen
Tel. +49 7461 96 48 55 - 0
Fax: +49 7461 96 48 55 - 95
Email: info@micromed.com
Web: www.micromed.com

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Tel +49 7461 96 48 550
Fax + 49 7461 96 48 55 95

info@micromed.com
www.micromed.com

© Copyright by Micromed Medizintechnik GmbH - Technical changes reserved - 06.2018