

Your advantage

- » Sharp, atraumatic and smooth
- » The fixation mandrin avoids unwanted tissue infections
- » The special cutting point and specialized surface treatment provide easy puncture even into deep tissue.
- » For better positioning the needles are marked with a galvanically applied not removable depth gauge, at intervals of 1 cm
- » Needle for biopsy and puncture for pain and tumor therapy

The logo for MICROMED features a stylized red wave above the brand name. The name 'MICROMED' is written in a red, spaced-out, sans-serif font. A horizontal line is positioned above the letters 'I', 'C', 'R', 'O', and 'M'.

M I C R O M E D[®]



Available: Diameter 18G and 23G
Lengths: 100 mm and 150 mm
Made in Germany

8,20 €
per piece

50 items per unit

CANNULAS
for pain and tumor therapy



Order sheet

Telephone

0049 7461 - 96 48 55-25

Fax

0049 7461 - 96 48 55-95

Mail

kfranco@micromed.com

Address: _____

FAO: _____

Street: _____

Postcode: City: _____

» In addition to surface polishing the cannulas/interventional micro probes (IMS) are dry siliconized so as to be able to reach deep lying tissue more easily

» The use of the special stainless steel mandrin enables precise punctuation to deeper tissue regions. This material is extremely strong and is combined with a slightly increased wall strength. The additionally siliconized surface completes the optimal needle

» Special, flexible stainless steel for easier correction of position, granting high stability

Description	Art.-No.	Diameter		Length	Package*
		[G]	ca. [mm]	[mm]	
Cannula 23 G with depth marker (blue)	EMS 23-10	23	0,65	100	
Cannula 18 G with depth marker (pink)	EMS 18-15	18	1,25	150	

*one pack contains 50 cannulas, price 410,- Euro, plus freight

Customer Service

Mrs. Franco 07461-96 48 55-25

Technical Service

Mr. Steidle 07461-96 48 55-21

CANNULAS
for pain and tumor therapy

