


125-501-018 180 mm
Absaug Instrument, gerade, biegsam
Suction-Coagulator, straight, malleable



125-502-018 180 mm
Absaug Instrument, gewinkelt
Suction-Coagulator, angled



125-503-011 110 mm
Absaug Instrument, 2x gewinkelt
Suction-Coagulator, 2x angled

Non-Stick Arbeitsende verhindert das Ankleben von Gewebe
Non-Stick tip to eliminate tissue sticking

Sicherheitsventil verhindert ein Festsaugen
Safety valves to prevent adherence of tissue

Kontrolle der Saugleistung
Suction control

Standard 4mm Anschluß
Standard 4mm electrical connector

Universeller Anschluß für Absaugschläuche
Universal connector for suction tubing